

# Need a Coughing Break? We are here to help!



## Cough Pause™ Option

The Vest® Systems, Model 105 and Model 205 devices have a **NEW** feature.

Our new Cough Pause™ Option allows you to program time into your treatment as a reminder to cough. The interval and duration options are programmable.



# ¿Habla Español?



## Languages

The Vest® Systems, Model 105 and Model 205 devices can be programmed in other languages. The device is programmed in English and one additional language can be added.

## Languages Available

- English
- German
- Italian
- French
- International Spanish
- Dutch
- Swedish
- International Portuguese